


Chemetron® Nitrogen Regulator Assembly

Product Specifications

- Lightweight
- Portable
- High flow (16-18 scfm)
- Wide operating range (5-180 psig)
- Variety of configurations

Performance Data

Operating range:	Adjustable 5-180 psig
Capacity:	Exceeds 16 scfm at 160 psig Outlet setting
	Exceeds 18 scfm at 180 psig Outlet setting



Description

Nitrogen regulator assemblies are designed to be used with the nitrogen outlet stations. Together they provide pneumatic instruments with varying pressures as required for operation.

The assembly includes an oversize adjustment knob, male Dyna-Con® or other fitting inlet adapter, 0-300 psig pressure gauge indicating the current regulated pressure, and an outlet valve.

Ordering Information

Catalog Number	Inlet	Outlet
11-20-0002	Dyna-Con Male, 90 Degree	Schrader Female Coupler (as shown above)
11-20-0003	Dyna-Con Male, 180 Degree	Schrader Female Coupler
11-20-0007	DISS Female, 180 Degree	Schrader Female Coupler
11-20-0021	DISS Female, 180 Degree	DISS Male w/Check
11-20-0022	DISS Female, 90 Degree	DISS Male w/Check
11-20-0023	Schrader Male, 180 Degree	Schrader Female Coupler
11-20-9045	DISS Female, 90 Degree	Schrader Female Coupler
11-20-9046	DISS Female, 180 Degree	1/4" NPT Male
11-20-9049	1/4" NPT Male, 180 Degree	1/4" NPT Female

Operating Instructions




1. Attach the regulator assembly to the turbo-pneumatic instrument.
2. Turn the regulator adjustment knob counter-clockwise to zero stop.
3. Plug the regulator into the nitrogen outlet and lock in place. Activate the flow from the outlet.
4. Adjust the regulator to the desired setting by turning the adjustment knob clockwise. Operate the tool. Adjust the regulator if necessary per the tool manufacturer's instructions.
5. When not in use, deactivate the flow at the Dyna-Con valve to bleed pressure from the regulator and hose. Then disassemble hose from the regulator if necessary.

Maintenance and Repair

1. Audibly inspect for leakage around the threaded connections before use. Do not use if leakage is detected.
2. Repair to the regulator itself is not economically feasible. See the Replacement Parts List for components.
3. Reassembly of threaded connections should include accepted sealing methods.

CAUTION: Always disconnect the nitrogen regulator from the outlet station before removing equipment from the regulator.

Replacement Parts



1.	11-91-0073	Dyna-Con Adapter x ¼ NPT
2.	12-90-0003	90° Elbow ¼ NPT X ¼ FPT
3.	11-20-9000	Regulator
4.	86-90-3009	0-300 PSI Gauge
5.	11-10-0095	Schrader Nitrogen Coupler
6.	11-20-0008	Dyna-Con check unit with 1/4 NPT Male (not illustrated)



Allied Healthcare Products, Inc. ▪ 1720 Sublette Avenue ▪ St. Louis, MO 63110 USA
314-771-2400 ▪ 800-444-3954 ▪ FAX: 314-771-0650 ▪ www.alliedhpi.com

U.S.
Toll Free: 800-444-3940
FAX: 800-477-7701

Canada & Mexico
Toll Free: 800-446-0552
FAX: 800-246-6201

International
Telephone: 314-268-1683
FAX: 314-771-5183



Allied Healthcare Products, Inc.—Global Support of Life

All specifications are nominal and
subject to change without notice.

Form 11-00-0003
Rev. MAY2016