

Multilator™



The LSP Multilator is a compact oxygen distribution system featuring eight LSP Constant Flow Selector Valves for precise delivery of oxygen to as many as nine patients from one cylinder.

The constant flow selector valves offer accurate and repeatable flow rates, regardless of position or visibility. All knobs are clearly marked and feature detents at each flow setting for accurate adjustment, even in low light conditions. The system comes complete with disposable masks and barbed fittings, and also can be ordered with optional bubble humidifiers. The custom case provides protection as well as storage space for accessories.

The Multilator is part of a sophisticated resuscitation/inhalation system that can expand as your needs grow.

Standard Features

- Easy to use ... no assembly needed.
- Durable, compact, and lightweight.
- Versatile ... provides administration of low flow oxygen and/or oxygen powered resuscitation.

Multilator™

SPECIFICATIONS

Inlet Pressure: 50 psig nominal

Operating Temperature: -30° F to 125° F

Storage Temperature: -40° F to 160° F

Outlets: 9 outlets, including 8 LSP Constant Flow Selector Valves with flow rates of 2, 4, 6, 10 and 15 lpm, and one standard male oxygen DISS check valve fitting

Filter: 25 micron stainless steel mesh screen on each selector valve

Masks: 8 adult size single-use masks, each with 7' of oxygen supply tubing

Hose Length: 10'

Weight: 12 lbs. (including humidifiers)

Material: Body - anodized aluminum; fittings - plated brass; case - h.d. polyethylene

ORDERING INFORMATION

L405 Multilator: 8 constant flow selector valves, 8 adult masks & tubing, 8 adapters, 1 10' oxygen hose, 8 humidifiers, case

L405-010 Multilator: 8 constant flow selector valves, 8 adult masks & tubing, 8 adapters, 1 10' oxygen hose, case

L233 Constant Flow Selector Valve, 1/8" female NPT inlet, oxygen DISS/barbed outlet, flow rates 2, 4, 6, 10, 15 lpm

L233-020 Constant Flow Selector Valve, 1/8" female inlet, barbed outlet, flow rates 2, 4, 6, 10, 15 lpm

L233-050 Constant Flow Selector Valve, 1/8" female inlet, oxygen DISS/barbed outlet, flow rates 2, 4, 6, 10, 26 lpm

L040114 Case for Multilator

64041 Oxygen Mask, adult, with 7' tubing, Latex free, 50/case

64375 Disposable Humidifier, 300 cc working volume, 3 psi, plastic nut, 50/case

65401 Oxygen Nut & Stem Adapter, 50/case



L233 Constant Flow Selector Valve shown with 65401 Nut & Stem Adapter



64375 Disposable Humidifier



64041 Oxygen Mask



EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands



ISO 13485: 2003

Allied Healthcare Products, Inc. ▪ 1720 Sublette Avenue ▪ St. Louis, MO 63110 USA
314/771-2400 ▪ 800/444-3954 ▪ FAX: 314/771-0650 ▪ www.alliedhpi.com

U.S.
Toll Free: 800/444-3940
FAX: 800/477-7701

Canada & Mexico
Toll Free: 800/446-0552
FAX: 800/246-6201

International
Telephone: 314/268-1683
FAX: 314/771-5183

©Allied Healthcare Products, Inc.®

All specifications are nominal and subject to change without notice.

Form 90-00-0015
Rev. May 2017