

# *Infant / Pediatric Immobilization Board*



MRI Compatible  
and  
X-Ray Translucent

The LSP Infant/Pediatric Immobilization Board is a versatile system capable of easily immobilizing infants and children with spinal injuries, or it can be used as a restraint device. It will accommodate children from infant to approximately 75 pounds (34 kg).

The unique harness system, which is integral to the device, allows easy manipulation of the child for proper airway management or radiological procedures. Once the child is positioned on the immobilization board, the spine is functionally maintained in a neutrally-aligned static position. The split-leg strapping design allows each leg to be immobilized and treated separately.

The removable and washable nylon covering is easy to clean and disinfect, and durable enough to withstand harsh treatment either in the field or in a hospital. The 1/2-inch board stiffener is made of finished Baltic birch, and is also easy to clean and maintain.

The Infant/Pediatric Board is completely radiotransparent and fully MRI compatible, allowing immobilization during x-ray or MRI.

The optional Drug Dosage and Patient Management Sleeve is designed to attach a Broselow Tape® or a customized patient management system as a reference for quicker determination of treatment.

# Infant / Pediatric Immobilization Board



L484C shown with optional  
Drug Dosage & Patient  
Management Sleeve (L484-075)



L720061-048 Carry Bag

## SPECIFICATIONS

**Length:** 47.5"

**Width:** 9.75"

**Depth:** 0.5"

**Weight:** 6 lbs.

**Capacity:** 75 lbs.

**Outer Covering:** 500 denier Cordura

**Backboard:** 0.5" Baltic birch

**Closures:** hook & loop

## ORDERING INFORMATION

|                    |  |
|--------------------|--|
| <b>L484C</b>       | Infant/Pediatric Immobilization Board with plastic clips               |
| <b>L484-040C</b>   | Infant/Pediatric Immobilization Board with plastic clips and carry bag |
| <b>L484-075</b>    | Drug Dosage and Patient Management Sleeve for L484C and L484-040C      |
| <b>L720061-048</b> | Carry Bag for L484C and L484-040C                                      |
| <b>L720073</b>     | Replacement Head Harness   |
| <b>L018077</b>     | Replacement Board Stiffener  |



L484-C  
L484-040C



**EMERGO EUROPE**  
Prinsessegracht 20  
2514 AP The Hague  
The Netherlands



Allied Healthcare Products, Inc. ▪ 1720 Sublette Avenue ▪ St. Louis, MO 63110 USA  
314/771-2400 ▪ 800/444-3954 ▪ FAX: 314/771-0650 ▪ [www.alliedhpi.com](http://www.alliedhpi.com)

U.S.  
Toll Free: 800/444-3940  
FAX: 800/477-7701

Canada & Mexico  
Toll Free: 800/446-0552  
FAX: 800/246-6201

International  
Telephone: 314/268-1683  
FAX: 314/771-5183

**ISO 13485: 2003**

©Allied Healthcare Products, Inc.®

All specifications are nominal and  
subject to change without notice.

Form 90-00-0016  
Rev. May 2017