

# CPR Mouth to Mask Resuscitator



L64186 Mouth to Mask Resuscitator

**LSP Mouth to Mask Resuscitators** are designed to give any rescuer a safe and effective method of providing ventilation. This single patient use device incorporates a high-flow one-way valve and a 3M Filtrete™ air filter with 99.8% bacteria filtration efficiency and 99.6% viral filtration efficiency to help protect the rescuer from bacteria or contaminants from the patient.

## Key Features

- Filter strategically located above the one-way valve to prevent patient secretions from saturating the filter
- Designed for use in EMS, hospital, home care and industrial markets
- Standard mask size accommodates adults and children
- Operation is not affected by extreme temperatures
- Configurations include a hard case for compact storage, oxygen inlet, head strap, gloves and wipes

## ORDERING INFORMATION

<b>L62183</b>	Face Mask with one-way valve and filter, 12/case
<b>L64184</b>	Face Mask with one-way valve and filter in hard white case, 12/case
<b>L64184X</b>	Face Mask with one-way valve and filter, gloves & wipes in hard white case, 12/case
<b>L64185</b>	Face Mask with one-way valve and filter, O <sub>2</sub> inlet & head strap, 12/case
<b>L64186</b>	Face Mask with one-way valve and filter, O <sub>2</sub> inlet & head strap in hard white case, 12/case
<b>L62050</b>	One-way valve with filter, 50/case



**EMERGO EUROPE**  
Prinsessegracht 20  
2514 AP The Hague  
The Netherlands

L64184  
L64186



# Flexible CPR Face Shield



L63002 Face Shield  
Resuscitator

**LSP Flexible CPR Face Shields** provide an inexpensive barrier offering protection for the rescuer while administering rescue breathing. This single patient use device incorporates a high-flow one-way valve and a 3M Filtrete™ air filter with 99.8% bacteria filtration efficiency and 99.6% viral filtration efficiency.

## Key Features

- Filter strategically located above the one-way valve to prevent patient secretions from saturating the filter
- Directions for use are printed onto the face shield
- Designed for EMS, hospital, home care and industrial markets for both adult and child
- Operation is not affected by extreme temperature

## ORDERING INFORMATION

**L63002** Flexible CPR Face Shield with one-way valve and filter, 50/case



Allied Healthcare Products, Inc. ▪ 1720 Sublette Avenue ▪ St. Louis, MO 63110 USA  
314/771-2400 ▪ 800/444-3954 ▪ FAX: 314/771-0650 ▪ [www.alliedhpi.com](http://www.alliedhpi.com)

U.S.  
Toll Free: 800/444-3940  
FAX: 800/477-7701

Canada & Mexico  
Toll Free: 800/446-0552  
FAX: 800/246-6201

International  
Telephone: 314/268-1683  
FAX: 314/771-5183

ISO 13485: 2003

©Allied Healthcare Products, Inc.®

All specifications are nominal and  
subject to change without notice.

Form 90-00-0030  
Rev. May 2017