


Replacement Hose Assemblies for Ceiling Columns and Hose Reels


Stainless Steel Manual Retractable Ceiling Column



DISS Female x DISS Female Replacement Hose Assembly

Catalog Number	Service
83-90-0110	Oxygen
83-90-0111	Vacuum
83-90-0112	Nitrous Oxide
83-90-0113	Medical Air
83-90-0114	Nitrogen
83-90-0115	WAGD
83-90-0116	Carbon Dioxide
83-90-0117	O ₂ /CO ₂ < 7%

Stainless Steel Vacuum Powered Retractable Ceiling Column



DISS Female x DISS Female Coiled Synflex Tubing Replacement Hose Assembly

Catalog #	Service
83-90-0001	Oxygen
83-90-0002	Vacuum
83-90-0003	Nitrous Oxide
83-90-0004	Medical Air
83-90-0005	Nitrogen
83-90-0006	Carbon Dioxide


Ordering Information

The following constitute 15-foot Replacement Hose Assemblies for Chemetron Hose Reel Assemblies. The Hose Assembly is complete with gooseneck DISS female nut and nipple, Chemetron Female Quick Connect Coupler, and Rubber Ball Stop.

Catalog Number	Service
82-90-0236	Oxygen
82-90-0237	Vacuum
82-90-0238	Nitrous Oxide
82-90-0239	Medical Air
82-90-0240	Nitrogen
82-90-0241	WAGD
82-90-0242	CO2

The following constitute Replacement Hose Assemblies for the supply hose (short length between the hose reel inlet and DISS outlet station) inside the Hose Reel Box. The Hose Assembly is complete with end fittings.

Catalog Number	Description
82-90-0225	Oxygen
82-90-0226	Vacuum
82-90-0227	Nitrous Oxide
82-90-0228	Medical Air
82-90-0229	Nitrogen
82-90-0230	WAGD
82-90-0232	CO2



NOTE: For any replacement parts for old style (pre 1977) 15-foot hose reel assemblies, contact the Service Department at 314-268-1600 or 800-411-5136.

All specifications are nominal and subject to change without notice.

Warranty
See Allied Statement of Warranties for details

© 2006 Allied Healthcare Product, Inc.® Printed in USA

Customer Service:
USA: (800) 444-3960
FAX: (314) 771-7537
Canada/Mexico: (800) 446-0552
FAX: (800) 246-6201
International: (314) 268-1683
FAX: (314) 771-5183



1720 Sublette Avenue
St. Louis, MO 63110 USA
(314)771-2400
www.alliedhpi.com
ISO 9001 CERTIFIED