

Gomco® 305 Aspirator

The features are new, but the Gomco quality is the same.



Model 3005, with 1500 ml disposable canister

The Gomco AC Model 305 complies with UL 60601-1 standards for use in hospitals, nursing home facilities and home care, and is licensed for use in Canada. This durable unit meets the need for controlled vacuum regulation up to 22" Hg with a flow rate of greater than 30 lpm. The aspirator has a high-impact plastic housing and dependable diaphragm pump, reduced operating noise level, and weighs less than 10 pounds. Each unit comes with a 1500 ml disposable canister, a suction tubing kit, and three hydrophobic/bacteria filters.

Standard Features

- Durable high-impact plastic case
- Vibration-free easy-to-read regulator gauge
- Wire basket accommodates a variety of canisters
- Quick-connect tubing for all connections
- Hydrophobic/bacteria filter eliminates pump flooding and prevents contaminant damage
- 3-year warranty

Gomco® 305 Aspirator



01-22-3005

SPECIFICATIONS

Vacuum Range:	Up to 22" Hg
Flow Rate:	Greater than 30 lpm
Sound Level:	Less than 60 dBA
Operating Temperature:	50° F to 104° F
Dimensions:	9.4" H x 7.5" D x 16.8" W (239 mm x 191 mm x 427 ml)
Weight:	9.8 lbs. (4.4 kg)
Electrical Requirements:	115 volt AC/60 Hz/2.0 amps maximum
UL Approval:	Tested to UL 60601-1
Warranty:	3 years

ORDERING INFORMATION

01-22-3005	Aspirator, 1500 ml canister, suction tubing kit, 3 filters, 115/60
01-22-0305	Aspirator, 1200 ml glass bottle, suction tubing kit, 3 filters, 115/60

ACCESSORIES

S1160BA-RPL	Disposable Canister, 800 cc (5/case)
20-08-0016	Disposable Canister, 1500 ml (16/case)
01-90-3497	Glass Bottle, 1200 ml (6/case)
S610100	Suction Tubing Kit (1 each 72" & 15" blue-tipped PVC tubing, 1/4" ID)
01-90-3928	Disposable Hydrophobic/Bacteria Filter, 1/4" barb x 1/4" barb (3/pkg.)



ISO 13485: 2003

Allied Healthcare Products, Inc. • 1720 Sublette Avenue • St. Louis, Missouri 63110
314/771-2400 • 800/444-3954 • Fax: 314/771-0650 • www.alliedhpi.com
International offices in Europe, Asia, and South America

U.S.
Toll Free: 800/444-3940
Fax: 800/477-7701

Canada & Mexico
Toll Free: 800/446-0552
Fax: 800/246-6201

International
Telephone: 314/268-1683
Fax: 314/771-5183