



Tabletop Aspirator Model 300



01-22-0300

The Gomco 300 is a heavy duty general use tabletop aspirator for hospitals, clinics and physician's offices. It generates up to -559 mm Hg of vacuum pressure and a 29 lpm flow rate at open flow.

Key Features

- Maintenance-free and oil-free diaphragm pump ensures maximum reliability
- Similar performance and quality of a mobile aspirator in a tabletop design
- Easy to store, the 300 takes up less than one square foot of space
- Disposable hydrophobic bacteria filter helps prevent fluid and aerosol contamination

Specifications

Vacuum Range: -25 to -559 mm Hg
 Flow Rate: 29 lpm at open flow
 Sound Level: Less than 60 dBA
 Pump Type: Oil-less diaphragm
 Electrical Cord: 6' hospital grade
 Dimensions: 13" W x 9" D x 12" H
 (30 x 23 x 30 cm)
 Weight: 19 lbs. (9 kg)
 UL/ETL Listed: No
 Warranty: 1 year
 Regulatory Approvals: FDA, HC

Model	Voltage	Glass Bottle	Disposable Canister
01-22-0300	115V/60 Hz	1 - 600 ml	
01-22-3001	115V/60 Hz		1 - 1100 ml
01-23-0300	220V/50 Hz	1 - 600 ml	

All models come with disposable suction tubing and 3 disposable hydrophobic bacteria filters.

01-90-2000	Disposable Suction Tubing Kit, one 18" & one 72" per pkg.
01-90-2454	Glass Collection Bottle, with rubber cap, short & long tube assembly, 600 ml
01-90-3100	Disposable Hydrophobic Bacteria Filter, 1/8" NPT x 1/4" hose barb (3/pkg.)
01-90-3695	Disposable Collection Canister, stem inlet & hydrophobic filter, 1100 ml (12/case)



MADE IN USA



Allied Healthcare Products, Inc. ▪ 1720 Sublette Avenue ▪ St. Louis, MO 63110 USA
 314/771-2400 ▪ 800/444-3954 ▪ FAX: 314/771-0650 ▪ www.alliedhpi.com

U.S.
 Toll Free: 800/444-3940
 FAX: 800/477-7701

Canada & Mexico
 Toll Free: 800/446-0552
 FAX: 800/246-6201

International
 Telephone: 314/268-1683
 FAX: 314/771-5183

ISO 13485: 2003

©Allied Healthcare Products, Inc.®

All specifications are nominal and subject to change without notice.

Form 21-00-0300
 July 2013