



General Use Aspirators

Models 3020, 3040



01-62-3020



01-62-3040

The Gomco 3020 and 3040 are general use aspirators. Each model generates up to -559 mm Hg of vacuum pressure and a 26 lpm flow rate at open flow. The Gomco 3020 comes with a 600 ml or an 1100 ml collection unit and is ideal for clinics, whereas the Gomco 3040 has a high capacity 2200 ml or 2800 ml collection unit and is better suited for hospitals.

Key Features

- Maintenance-free and oil-free diaphragm pump ensures maximum reliability
- Easy-to-store and transport compact design with a flat stainless steel top
- Convenient storage for a 6' hospital grade electrical cord
- Disposable hydrophobic bacteria filter helps prevent fluid and aerosol contamination

Specifications

Vacuum Range: -25 to -559 mm Hg
 Flow Rate: 26 lpm at open flow
 Sound Level: Less than 52 dBA
 Pump Type: Oil-less diaphragm
 Electrical Cord: 6' hospital grade
 Dimensions: 18" W x 17" D x 33" H
 (46 x 43 x 84 cm)
 Weight: 52 lbs. (24 kg)
 UL/ETL Listed: Yes
 Warranty: 1 year
 Regulatory Approvals: FDA, HC

| Model | Voltage | Glass Bottle | Poly Bottle | Disposable Canister |
|------------|------------|--------------|-------------|---------------------|
| 01-62-3020 | 115V/60 Hz | 1 - 600 ml | | |
| 01-63-3020 | 220V/50 Hz | 1 - 600 ml | | |
| 01-62-3021 | 115V/60 Hz | | | 1 - 1100 ml |
| 01-62-3040 | 115V/60 Hz | | 1 - 2800 ml | |
| 01-62-3041 | 115V/60 Hz | | | 1 - 2100 ml |
| 01-63-3040 | 220V/50 Hz | | 1 - 2800 ml | |
| 01-64-3040 | 220V/60 Hz | | 1 - 2800 ml | |

All models come with disposable suction tubing and 3 disposable hydrophobic bacteria filters.

| | |
|------------|---|
| 01-90-2000 | Disposable Suction Tubing Kit, one 18" & one 72" per pkg. |
| 01-90-2454 | Glass Collection Bottle, with rubber cap, short & long tube assembly, 600 ml |
| 01-90-3100 | Disposable Hydrophobic Bacteria Filter, 1/8" NPT x 1/4" hose barb (3/pkg.) |
| 01-90-3694 | Polycarbonate Collection Bottle, with metal cover and cap & float assembly, 2800 ml |
| 01-90-3695 | Disposable Collection Canister, stem inlet and hydrophobic bacteria filter, 1100 ml (12/case) |
| 01-90-3712 | Disposable Collection Canister, stem inlet and hydrophobic bacteria filter. 2100 ml (10/case) |



MADE IN USA



Allied Healthcare Products, Inc. • 1720 Sublette Avenue • St. Louis, MO 63110 USA
 314/771-2400 • 800/444-3954 • FAX: 314/771-0650 • www.alliedhpi.com

U.S.
 Toll Free: 800/444-3940
 FAX: 800/477-7701

Canada & Mexico
 Toll Free: 800/446-0552
 FAX: 800/246-6201

International
 Telephone: 314/268-1683
 FAX: 314/771-5183

ISO 13485: 2003

©Allied Healthcare Products, Inc.®

All specifications are nominal and subject to change without notice.

Form 21-00-3020-3040
 July 2013