



# Surgical Aspirators

## Models 3910, 3940



01-22-3910



01-22-3940

Gomco 3910 and 3940 are heavy duty, high performance mobile aspirators engineered for surgical suites. Both models generate up to -635 mm Hg of vacuum pressure and a 34 lpm flow rate at open flow. The Gomco 3910 comes with one reusable collection bottle, and the Gomco 3940 comes with two selectable collection bottles for extended use.

### Key Features

- Maximum reliability thanks to the robust, maintenance-free and oil-free diaphragm pump design
- Stable, heavy duty cart design with flat stainless steel sides and top
- Mechanical shut-off float helps prevent fluid contamination
- Convenient storage drawer for suction accessories

### Specifications

Vacuum Range: -25 to -635 mm Hg  
 Flow Rate: 34 lpm at open flow  
 Sound Level: Less than 54 dBA  
 Pump Type: Oil-less belt-driven diaphragm  
 Electrical Cord: 6' hospital grade  
 Dimensions: 20" W x 15" D x 32" H (50 x 38 x 80 cm)  
 Weight: 110 lbs. (50 kg)  
 UL/ETL Listed: No  
 Warranty: 1 year  
 Regulatory Approvals: FDA

Model	Voltage	Poly Bottle
01-22-3910	115V/50 Hz	1 - 2800 ml
01-22-3940	115V/50 Hz	2 - 2800 ml
01-23-3940	220V/50 Hz	2 - 2800 ml

**Explosion-proof models 3810 and 3840 are also available. For information, contact your Allied Healthcare sales representative or Customer Service.**

*All models come with disposable suction tubing and 3 disposable hydrophobic bacteria filters.*

01-90-2000	Disposable Suction Tubing Kit, one 18" & one 72" per pkg.
01-90-3694	Polycarbonate Collection Bottle, with metal cover and cap & float assembly, 2800 ml



MADE IN USA



Allied Healthcare Products, Inc. ▪ 1720 Sublette Avenue ▪ St. Louis, MO 63110 USA  
 314/771-2400 ▪ 800/444-3954 ▪ FAX: 314/771-0650 ▪ www.alliedhpi.com

U.S.  
 Toll Free: 800/444-3940  
 FAX: 800/477-7701

Canada & Mexico  
 Toll Free: 800/446-0552  
 FAX: 800/246-6201

International  
 Telephone: 314/268-1683  
 FAX: 314/771-5183

ISO 13485: 2003

©Allied Healthcare Products, Inc.®

All specifications are nominal and subject to change without notice.

Form 21-00-3910-3940  
 July 2013