

Uterine Aspirators Models 3930, 4072



The Gomco 3930 and 4072 are mobile aspirators specialized for dilation and curettage (D&C). Both models generate up to -635 mm Hg of vacuum pressure. The Gomco 3930 is our heavy duty hospital model, and the Gomco 4072 is our clinic model.

Key Features

- Two 600 ml glass collection bottles with a specimen collection bag
- Comes with four reusable curettes (unless otherwise noted)
- Chemical-resistant easy to clean stainless steel top
- Convenient storage for cord with hospital-grade electrical plug
- Disposable hydrophobic bacterial filter helps prevent fluid and aerosol contamination (Models 4072/4073)

Specifications

Vacuum Range:

Model 3930: -25 to -635 mm Hg Model 4072: -80 to -635 mm Hg

Flow Rate:

Model 3930: 34 lpm at open flow Model 4072: 40 lpm at open flow Sound Level: Less than 54 dBA

Pump Type:

Model 3930: Oil-less belt-driven diaphragm

Model 4072: Carbon rotary vane Electrical Cord: 6' hospital grade

Dimensions:

Model 3930: 20" W x 15" D x 32" H

(50 x 38 x 80 cm)

Model 4072: 18" W x 17" D x 33" H

(46 x 43 x 84 cm)

Weight:

Model 3930: 110 lbs. (50 kg) Model 4072: 46 lbs. (21 kg)

UL/ETL Listed:

Model 3930: No Model 4072: Yes Warranty: 1 year

Regulatory

Approvals: FDA, HC (model 4072 only)

Model	Voltage	Glass Bottle	Curettes
01-22-3930	115V/50-60 Hz	2 - 600 ml	4 - 1 of each size
01-52-4072	115V/50-60 Hz	2 - 600 ml	4 - 1 of each size
01-52-4073	115V/50-60 Hz	2 - 600 ml	

Explosion-proof model 3830 is also available. For information, contact your Allied Healthcare sales representative or Customer Service.

Model 3930 comes with disposable suction tubing (01-90-2235) and no bacteria filters. Models 4072/4073 come with disposable suction tubing (01-90-2000) and 3 disposable hydrophobic bacteria filters.

01-90-2000	Disposable Suction Tubing Kit, one 18" & one 72" per pkg.		
01-90-2235	Disposable Suction Tubing, PVC, clear, 1/2" ID x 72"		
01-90-2675	Glass Collection Bottle (bottle only), 600 ml		
01-90-3010	Top, rubber, for 600 ml collection bottle		
01-90-3100	Disposable Hydrophobic Bacteria Filter, 1/8" NPT x 1/4" hose bard (3/pkg.)		
01-90-3078	Curette Handle		
01-90-8820	Curette Tip, 9.75" long, 6 mm OD		
01-90-8830	Curette Tip, 9.75" long, 9 mm OD		
01-90-8840	Curette Tip, 9.75" long, 11 mm OD		
01-90-8850	Curette Tip, 9.75" long, 13 mm OD		
01-90-9035	Disposable Suction Tubing, clear PVC, 6" x 1/4" ID		





Allied Healthcare Products, Inc. • 1720 Sublette Avenue • St. Louis, MO 63110 USA 314/771-2400 • 800/444-3954 • FAX: 314/771-0650 • www.alliedhpi.com

U.S. Toll Free: 800/444-3940 FAX: 800/477-7701 Canada & Mexico Toll Free: 800/446-0552 FAX: 800/246-6201 International Telephone: 314/268-1683 FAX: 314/771-5183

ISO 13485: 2003

All specifications are nominal and subject to change without notice.

Form 21-00-3930-4072 July 2013