

Schuco S5200 Pediatric Nebulizer



Introducing the Schuco® S5200 Pediatric Nebulizer ... Providing children the treatment they depend on with the comfort they deserve.

- Child-friendly design
- Lightweight yet quiet
- Easy to use and clean

Complete with all that's needed to begin the child's nebulizer therapy:

- Panda compressor
- Nebulizer cup
- Angled mouthpiece and mask
- Over six feet of oxygen tubing
- Five extra compressor filters
- 5-year limited warranty

**Schuco**
byAllied

Schuco® S5200

STANDARD FEATURES

- Quiet operation ... less than 60 dBA
- Lightweight ... just 3-1/2 lbs.
- Substantial compressor pressure ... 35-45 psig
- High compressor operating flow rate ... 7-11 lpm
- Therapeutic particle size ... 0.5 µm-5µm
- AC operation ... 120V 60 Hz
- UL classified
- 5-year limited warranty



Each S5200 is packaged in a retail-ready box.

ORDERING INFORMATION

S5200 Panda Compressor, nebulizer cup, angled mouthpiece, mask, 6 feet of oxygen tubing, 5 extra compressor filters

ACCESSORIES

S66-0126 Replacement Filters, 5/pkg. (20 pkgs./case)

64095 Pediatric Aerosol Mask with nebulizer and 7' clear tubing (50/case)

64085 Adult Aerosol Mask with nebulizer and 7' clear tubing (50/case)

61399 Replacement Nebulizer Kit: nebulizer cup, 7' clear tubing, mouthpiece, tee, 6" flex tube (50/case)

61400 Replacement Nebulizer Kit: nebulizer cup, 7' clear tubing, mouthpiece, tee (50/case)

88100 AC to DC Power Converter for use in car



Allied Healthcare Products, Inc. • 1720 Sublette Avenue • St. Louis, MO 63110 USA
314/771-2400 • 800/444-3954 • FAX: 314/771-0650 • www.alliedhpi.com

U.S.
Toll Free: 800/444-3940
FAX: 800/477-7701

Canada & Mexico
Toll Free: 800/446-0552
FAX: 800/246-6201

International
Telephone: 314/268-1683
FAX: 314/771-5183

ISO 13485: 2003

©Allied Healthcare Products, Inc.®

All specifications are nominal and subject to change without notice.

Form 23-00-0006
Rev. February 2011