

Chemetron® Medical Gas Outlet Repair Kit



The solution to your medical gas outlet maintenance needs is in this kit.

Your medical gas outlets must perform when you need them. That means maintaining them regularly. With the Chemetron® Medical Gas Outlet Repair Kit, your maintenance staff is equipped with the necessary repair parts to provide any needed outlet service, quickly and efficiently.

Most quick-connect medical gas outlets use an O-ring to seal around the adapter for a gas-tight connection. Occasionally that O-ring may need replacement. Chemetron offers a case for quick repair.

Supporting our case is the unique Chemetron O-ring service tool. It is the only tool of its kind in the industry that permits easy replacement of quick-connect outlet station O-rings – without disassembly of the outlet, without interrupting adjacent outlets, without shutting off the gas supply.

In the Chemetron Medical Gas Outlet Repair Kit, you will find a common parts kit with O-rings for all Chemetron outlet stations and components for 348, 378 and 400 series, quick-connect, and DISS outlets. Also contained in the complete repair kit is the innovative O-ring service tool and parts kit for quick-connect, DISS and high-pressure nitrogen outlet stations.

If you do not require all of these parts kits, ask about our starter kit.

The evidence is clear. There is no other medical gas outlet repair kit like this. Anywhere.

Let Chemetron solve your maintenance needs. Contact us today for more information or a demonstration.



Chemetron Medical Gas Outlet Repair Kit

Starter Repair Kit allows you to tailor the repair kit to your needs and budget. It is recommended that the starter kit be purchased for basic repair of outlet stations except high-pressure nitrogen.

Complete Repair Kit gives you all of the repair kits under one catalog number for convenience in ordering.

For complete repair of an outlet by specific type and model (quick-connect, DISS, and/or high-pressure nitrogen), order the respective parts kits as well. Space is provided in the starter kit carrying case for addition of the parts kits.

Order This Now	Catalog #	Provides You With	Allows You To Service or Repair
Starter Repair Kit	64-90-2023	Carrying case, labeled common parts kit, O-ring lubricant, O-ring service tool, and space for parts kits listed below	O-rings for all Chemetron (NCG®) quick-connect and DISS outlet stations, base housings, vacuum springs, and secondary checks for 348/346, 378/376, and 400/460 series outlets
Order Now or Later	Catalog #	Provides You With	Allows You To Service or Repair
Quick-connect 400 series outlet station parts kits	64-90-2001	Labeled organizer kit box, valve bodies, poppets, latch tips and buttons, latch springs	Series 348, 378, and 400 quick-connect outlet station valves and faceplates
DISS 460 series outlet station parts kit	64-90-2002	Labeled organizer kit box, valve bodies for most commonly used gas services, poppets, dust caps, gas service buttons and screws	Series 346, 376 and 460 DISS outlet station valves and faceplates
Quick-connect 500 series outlet station parts kits	64-90-2005	Labeled organizer kit box, secondary check poppets, tool, valve bodies, springs	Series 500 quick-connect outlet stations
DISS 560 series outlet station parts kit	64-90-2010	Labeled organizer kit box, secondary check poppets, tool, valve bodies, springs for most commonly used gas services	Series 560 DISS outlet stations
High-pressure nitrogen outlet station parts kit	64-90-2004	Labeled organizer kit box, Dyna-Con valve assemblies, housings, Dyna-Con service tool, knobs, lock-washers, O-rings and seals	Chemetron high-pressure Dyna-Con outlet station valves for nitrogen service

Order This Now	Catalog #	Provides You With	Allows You To Service or Repair
Complete Repair Kit	64-90-2028	Starter Repair Kit Quick-Connect, DISS and High-Pressure Nitrogen Parts Kit	All of the outlet stations listed above except 500 and 560 series outlets

All specifications are nominal and subject to change without notice.

Warranty:
See Allied Statement of Warranties for details.

© 2005 Allied Healthcare Products, Inc.® Printed in USA

Customer Service:

USA: (800) 444-3960
 FAX: (314) 771-7537
 Canada/Mexico: (800) 446-0552
 FAX: (800) 246-6201
 International: (314) 268-1683
 FAX: (314) 771-5183



1720 Sublette Avenue
 St. Louis, MO 63110 USA
 (314)771-2400
 www.alliedhpi.com
 ISO 9001 CERTIFIED

Product Specification