


The HDx Backboard®

The HDx is a leading composite spine immobilization board known for its unsurpassed rigidity and lightweight construction.

The seamless HDPE shell is fluid-resistant, and the composite internal stiffening system provides unmatched rigidity with minimal x-ray obstruction.

- Strong and durable, yet extremely lightweight at only 13.6 lbs
- Low profile design is only 1.5" thick
- 10 large elevated hand-holds for easy patient handling
- Separate strap holes for easy immobilization
- One-piece, seamless HDPE shell is impervious to blood and other substances for high resistance to contamination, corrosion or decay



SPECIFICATIONS

Length:	72" (182.88 cm)
Width:	16" (40.64 cm)
Depth:	1.5" (3.81 cm) including runners
Weight:	13.6 lbs. (6 kg)
Capacity:	500 lbs. (226 kg)
Material:	High-density polymer

ORDERING INFORMATION

L712000	HDx, yellow, no pins
L712010	HDx, yellow, 10 pins
L715000	HDx, red, no pins
L715010	HDx, red, 10 pins
L716000	HDx, blue, no pins
L716010	HDx, blue, 10 pins
L718000	HDx, orange, no pins
L718010	HDx, orange, 10 pins

Can be imprinted with large-format custom graphics (22.5" wide x 7.5" high area) with no setup charges. Contact Allied's Customer Service Department for assistance with ordering custom graphics.



MADE IN USA



Allied Healthcare Products, Inc. ▪ 1720 Sublette Avenue ▪ St. Louis, MO 63110 USA
314/771-2400 ▪ 800/444-3954 ▪ FAX: 314/771-0650 ▪ www.alliedhpi.com

U.S.
Toll Free: 800/444-3940
FAX: 800/477-7701

Canada & Mexico
Toll Free: 800/446-0552
FAX: 800/246-6201

International
Telephone: 314/268-1683
FAX: 314/771-5183

ISO 13485: 2003

©Allied Healthcare Products, Inc.®

All specifications are nominal and subject to change without notice.

Form 90-00-0013
Rev. March 2014