

Minilator™



L419-010

The LSP Minilator is a simple to operate constant flow device which provides oxygen for up to five patients simultaneously.

The versatile Minilator features 50 psig outlets, which can be used to operate a demand resuscitator valve, suction device, CPAP mask, or ventilator. All five oxygen outlets can function simultaneously, supplying 40 to 90 psig or low flow oxygen. Threaded barbed fittings with metered orifices can be fitted to the DISS check valves to provide consistent accurate flow rates of 4, 6, 10 or 15 liters per minute. Using standard DISS outlets, the Minilator is easily attached to an oxygen source via standard oxygen hoses and/or connectors. The zippered nylon carry case keeps the components clean and stores them fully assembled and ready for use.

Standard Features

- Expandable ... Minilators can be connected in tandem for emergency or mass casualty situations where more than five oxygen outlets are required.
- Durable ... anodized aluminum construction will withstand harsh weather conditions and corrosive materials.
- Compact and lightweight ... low-profile design for easy storage.

Minilator™

SPECIFICATIONS

Supply Pressure: 40-90 psig
Operating Temperature: -30° F to 125° F
Storage Temperature: -40° F to 160° F
Flow Rates: 4, 6, 10 or 15 lpm, depending on orifice
Materials: manifold - anodized aluminum;
fittings - plated brass;
barbs - H.I. polystyrene;
orifices - stainless steel
Weight: 13.5 oz.
Dimensions: 1" x 3" x 6" without barbs
1" x 6" x 7-1/2" with barbs attached

ORDERING INFORMATION

L419-010 Minilator: 4" aluminum manifold with five oxygen DISS check valves, five screw-on 4 lpm orifices, zippered nylon carry case
L419-020 Minilator, same as L419-010 but with five 6 lpm orifices
L419-030 Minilator, same as L419-010 but with five 10 lpm orifices
L419-050 Minilator, same as L419-010 but with five 15 lpm orifices
L419 Minilator, 4" manifold only
L555005 Orifice, 2 lpm, 5/pkg
L555006 Orifice, 4 lpm, 5/pkg.
L555007 Orifice, 6 lpm, 5/pkg.
L555008 Orifice, 10 lpm, 5/pkg.



L419-010



EMERGO EUROPE
Molenstraat 15
2513 BH, The Hague
The Netherlands



Allied Healthcare Products, Inc. • 1720 Sublette Avenue • St. Louis, MO 63110 USA
314/771-2400 • 800/444-3954 • FAX: 314/771-0650 • www.alliedhpi.com

U.S.
Toll Free: 800/444-3940
FAX: 800/477-7701

Canada & Mexico
Toll Free: 800/446-0552
FAX: 800/246-6201

International
Telephone: 314/268-1683
FAX: 314/771-5183

ISO 13485: 2003

©Allied Healthcare Products, Inc.®

All specifications are nominal and subject to change without notice.

Form 90-00-0014
Rev. January 2014