

# Multilator™



The LSP Multilator is a compact oxygen distribution system featuring eight LSP Constant Flow Selector Valves for precise delivery of oxygen to as many as nine patients from one cylinder.

The constant flow selector valves offer accurate and repeatable flow rates, regardless of position or visibility. All knobs are clearly marked and feature detents at each flow setting for accurate adjustment, even in low light conditions. The system comes complete with disposable masks and barbed fittings, and also can be ordered with optional bubble humidifiers. The custom case provides protection as well as storage space for accessories.

The Multilator is part of a sophisticated resuscitation/inhalation system that can expand as your needs grow.

## Standard Features

- Easy to use ... no assembly needed.
- Durable, compact, and lightweight.
- Versatile ... provides administration of low flow oxygen and/or oxygen powered resuscitation.

# Multilator™

## SPECIFICATIONS

**Inlet Pressure:** 50 psig nominal

**Operating Temperature:** -30° F to 125° F

**Storage Temperature:** -40° F to 160° F

**Outlets:** 9 outlets, including 8 LSP Constant Flow Selector Valves with flow rates of 2, 4, 6, 10 and 15 lpm, and one standard male oxygen DISS check valve fitting

**Filter:** 25 micron stainless steel mesh screen on each selector valve

**Masks:** 8 adult size single-use masks, each with 7' of oxygen supply tubing

**Hose Length:** 10'

**Weight:** 12 lbs. (including humidifiers)

**Material:** Body - anodized aluminum; fittings - plated brass; case - h.d. polyethylene

## ORDERING INFORMATION

**L405** Multilator: 8 constant flow selector valves, 8 adult masks & tubing, 8 adapters, 1 10' oxygen hose, 8 humidifiers, case

**L405-010** Multilator: 8 constant flow selector valves, 8 adult masks & tubing, 8 adapters, 1 10' oxygen hose, case

**L233** Constant Flow Selector Valve, 1/8" female NPT inlet, oxygen DISS/barbed outlet, flow rates 2, 4, 6, 10, 15 lpm

**L233-020** Constant Flow Selector Valve, 1/8" female inlet, barbed outlet, flow rates 2, 4, 6, 10, 15 lpm

**L233-050** Constant Flow Selector Valve, 1/8" female inlet, oxygen DISS/barbed outlet, flow rates 2, 4, 6, 10, 26 lpm

**L040114** Case for Multilator

**64041** Oxygen Mask, adult, with 7' tubing, Latex free, 50/case

**64375** Disposable Humidifier, 300 cc working volume, 3 psi, plastic nut, 50/case

**65401** Oxygen Nut & Stem Adapter, 50/case



L233 Constant Flow Selector Valve shown with 65401 Nut & Stem Adapter



64375 Disposable Humidifier



64041 Oxygen Mask



**EMERGO EUROPE**  
Molenstraat 15  
2513 BH, The Hague  
The Netherlands



Allied Healthcare Products, Inc. ▪ 1720 Sublette Avenue ▪ St. Louis, MO 63110 USA  
314/771-2400 ▪ 800/444-3954 ▪ FAX: 314/771-0650 ▪ www.alliedhpi.com

U.S.  
Toll Free: 800/444-3940  
FAX: 800/477-7701

Canada & Mexico  
Toll Free: 800/446-0552  
FAX: 800/246-6201

International  
Telephone: 314/268-1683  
FAX: 314/771-5183

ISO 13485: 2003

©Allied Healthcare Products, Inc.®

All specifications are nominal and subject to change without notice.

Form 90-00-0015  
Rev. January 2014