

Infant / Pediatric Immobilization Board



L484C shown with optional L484-075 Drug Dosage & Patient Management Sleeve (Broselow[®] Pediatric Emergency Tape not included)

MRI Compatible
and
X-Ray Translucent

The LSP Infant/Pediatric Immobilization Board is a versatile system capable of easily immobilizing infants and children with spinal injuries, or it can be used as a restraint device. It will accommodate children from infant to approximately 75 pounds (34 kg).

The unique harness system, which is integral to the device, allows easy manipulation of the child for proper airway management or radiological procedures. Once the child is positioned on the immobilization board, the spine is functionally maintained in a neutrally-aligned static position. The split-leg strapping design allows each leg to be immobilized and treated separately.

The removable and washable nylon covering is easy to clean and disinfect, and durable enough to withstand harsh treatment either in the field or in a hospital. The 1/2-inch board stiffener is made of finished Baltic birch, and is also easy to clean and maintain.

The Infant/Pediatric Board is completely radiotransparent and fully MRI compatible, allowing immobilization during x-ray or MRI.

The optional Drug Dosage and Patient Management Sleeve is designed to attach a Broselow Tape[®] or a customized patient management system as a reference for quicker determination of treatment.

Infant / Pediatric Immobilization Board



L484C shown with optional
Drug Dosage & Patient
Management Sleeve (L484-075)



L720061-048 Carry Bag

SPECIFICATIONS

Length: 47.5"

Width: 9.75"

Depth: 0.5"

Weight: 6 lbs.

Capacity: 75 lbs.

Outer Covering: 500 denier Cordura

Backboard: 0.5" Baltic birch

Closures: hook & loop

ORDERING INFORMATION

L484C	Infant/Pediatric Immobilization Board with plastic clips
L484-040C	Infant/Pediatric Immobilization Board with plastic clips and carry bag
L484-075	Drug Dosage and Patient Management Sleeve for L484C and L484-040C
L720061-048	Carry Bag for L484C and L484-040C
L720073	Replacement Head Harness
L018077	Replacement Board Stiffener



L484-C
L484-040C



EMERGO EUROPE
Molenstraat 15
2513 BH, The Hague
The Netherlands



Allied Healthcare Products, Inc. ▪ 1720 Sublette Avenue ▪ St. Louis, MO 63110 USA
314/771-2400 ▪ 800/444-3954 ▪ FAX: 314/771-0650 ▪ www.alliedhpi.com

U.S.
Toll Free: 800/444-3940
FAX: 800/477-7701

Canada & Mexico
Toll Free: 800/446-0552
FAX: 800/246-6201

International
Telephone: 314/268-1683
FAX: 314/771-5183

ISO 13485: 2003

©Allied Healthcare Products, Inc.®

All specifications are nominal and
subject to change without notice.

Form 90-00-0016
Rev. January 2014