

# *Miller Full Body Splint / Litter*



The Miller Full Body Splint/Litter is a versatile spine immobilization board and Stokes basket-compatible extrication litter. The system comes complete with a full head harness and integrated immobilization straps.

This unique harness system with integrated straps is ideal for water rescue as well as any rescue situation where non-integrated straps could become entangled or be lost. An optional head harness allows this same critical immobilization of a helmeted victim.

When used in conjunction with the LSP Half Back, the Miller board can be used for safe vertical lifting in normal and confined space rescue without compromising patient safety or spine immobilization.

The many unique features of the Miller board have made this the most popular rescue and immobilization product for more than 20 years.

- Rigid foam core provides both strength and buoyancy for water rescue immobilization
- Will float 250 pounds in fresh water and 300 pounds in salt water
- Integrated strapping system allows the rescuer to quickly and efficiently immobilize and restrain combative patients and avoids the loss of the straps
- Split-leg design allows for easy independent leg traction, accessibility, and compatibility with the Stokes litter basket
- Shoulder strap design prevents both lateral and horizontal movement of the patient during extrication and transportation

# Miller Full Body Splint/Litter



## SPECIFICATIONS

**Length:** 65.75" (167 cm)

**Width:** 14" (35.56 cm)

**Depth:** 3.25" (8.25 cm) at fullest point

**Weight:** 15 lbs. (6.81 kg)

**Load Capacity:** 500 lbs. (227 kg)

**Flotation Capacity:** 250 lbs. (114 kg) in fresh water; 300 lbs. (136 kg) in salt water

**Outer Plastic Shell:** High density polyethylene; temperature range -180°F to 240°F; resistant to gasoline and other chemicals; clean with soap and water

**Injected Foam Core:** Rigid polyurethane foam; meets military specifications.

**Webbing:** #7226 commercial-grade nylon, 2"; 5,000/6,000 lbs. (2700 kg) test; washable

**Closures:** 2" (5 cm) black hook and loop; washable; meets military specifications

**Nylon Fabric:** 100% utility nylon, black, 200 denier

**Foam:** #2 ethafoam; meets military specifications

**Slide:** Nickel-plated steel

## ORDERING INFORMATION

- L700** Miller Full Body Splint/Litter with V head rest & standard head harness  
**L701** Miller Full Body Splint/Litter with V head rest & head harness, military (olive)  
**L702** Miller Full Body Splint/Litter with V head rest & head harness, military (yellow)  
**L720001** Replacement V Head Rest for L700, L701 & L702 & L710  
**L720003-010** Helmet Immobilizer for L700, L701 & L702  
**L720016** Carry Case for L700, L701 & L702  
**L720077** Replacement Head Immobilizer for L700, L701 & L702  
**L720-015** 12" Strap Extension Kit for L700, L701 & L702 (6/kit)



L700  
L701  
L702



**EMERGO EUROPE**  
Molenstraat 15  
2513 BH, The Hague  
The Netherlands



Allied Healthcare Products, Inc. ▪ 1720 Sublette Avenue ▪ St. Louis, MO 63110 USA  
314/771-2400 ▪ 800/444-3954 ▪ FAX: 314/771-0650 ▪ [www.alliedhpi.com](http://www.alliedhpi.com)

U.S.  
Toll Free: 800/444-3940  
FAX: 800/477-7701

Canada & Mexico  
Toll Free: 800/446-0552  
FAX: 800/246-6201

International  
Telephone: 314/268-1683  
FAX: 314/771-5183

ISO 13485: 2003

©Allied Healthcare Products, Inc.®

All specifications are nominal and subject to change without notice.

Form 90-00-0024  
Rev. January 2014