

# *L238 Reusable Manual Resuscitator*



L238-220 Adult  
Reusable Manual  
Resuscitator

***L238 Reusable Manual Resuscitators*** are front-fill manual resuscitation systems that are available in adult and child sizes. They are fully autoclavable at 134°C (273°F) or cold sterilizable. The L238 features a silicone bag which provides an excellent “feel” of the patient’s respiratory response, and the durable polysulfone patient valve and bag adapter ensure an extended life even when subjected to repeated autoclave sterilization cycles.

The L238 comes with our patented Tru-Fit mask which molds to the patient’s face while maintaining a secure fit. This flexible and durable mask retains its original shape and is autoclavable at 134°C (273°F).

When connected to a metered low flow oxygen source via the L003482 oxygen enrichment assembly, the L238 resuscitator can deliver FiO<sub>2</sub> concentrations between 73% and 97%. Used in conjunction with an L471 bag refill valve, it ensures 100% oxygen is delivered during resuscitation.

The L238 resuscitator can be ordered in individual sizes or in combination kits that contain adult and child bags and Tru-Fit masks.

# Reusable Manual Resuscitator

<b>L238</b>	<b>Tru-Fit Mask Size</b>	<b>Oxygen Enrichment Assembly</b>	<b>Plastic Carrying Case</b>
<b>ADULT</b>			
L238-220*	Size 5		
L238-222*	Size 5	✓	
<b>CHILD</b>			
L238-210*	Size 2		
L238-212*	Size 2	✓	
<b>ADULT &amp; CHILD KIT</b>			
L238-300*	Sizes 5 & 2	✓	✓
<b>ACCESSORIES</b>			
LPEEP	Disposable PEEP Valve with adapter, 12/case		
L003482	Oxygen Enrichment Assembly		
L040081	Replacement orange molded plastic carry case (order separately)		
*CE 0413			

**Bag Volume:** Adult: 1700 ml  
Child: 680 ml  
Infant: 370 ml

**Patient Valve Dead Space:** 8 ml (not including mask)

**Patient Connection:** 22 mm outside diameter  
15 mm inside diameter

**PEEP Connection:** 30 mm outside diameter

**Operating Temperature:** -30°F to 125°F (-34°C to 52°C)

**Storage Temperature:** -40°F to 160°F (-40°C to 71°C)

**Sterilization Temperature:** 134°C

**Materials**

**Bag:** Silicone Rubber  
**Patient Valve:** Polysulfone  
**Bag Adapter:** Polysulfone  
**Mask:** Silicone Rubber



Optional disposable 30 mm PEEP valve with 19 mm valve adapter to complete the resuscitation system.

- Adjustable from 0-20 cm H<sub>2</sub>O
- Fits easily on the exhalation valve
- Can be used on any make of manual resuscitator



L238-222 Adult Reusable Manual Resuscitator with oxygen enrichment assembly



L040081 Replacement Carry Case



**EMERGO EUROPE**  
Molenstraat 15  
2513 BH, The Hague  
The Netherlands



Allied Healthcare Products, Inc. • 1720 Sublette Avenue • St. Louis, MO 63110 USA  
314/771-2400 • 800/444-3954 • FAX: 314/771-0650 • www.alliedhpi.com

U.S.  
Toll Free: 800/444-3940  
FAX: 800/477-7701

Canada & Mexico  
Toll Free: 800/446-0552  
FAX: 800/246-6201

International  
Telephone: 314/268-1683  
FAX: 314/771-5183

ISO 13485: 2003

©Allied Healthcare Products, Inc.®

All specifications are nominal and subject to change without notice.

Form 90-00-0029  
Rev. January 2014