

CPR Mouth to Mask Resuscitator



L64186 Mouth to Mask Resuscitator

LSP Mouth to Mask Resuscitators are designed to give any rescuer a safe and effective method of providing ventilation. This single patient use device incorporates a high-flow one-way valve and a 3M Filtrete™ air filter with 99.8% bacteria filtration efficiency and 99.6% viral filtration efficiency to help protect the rescuer from bacteria or contaminants from the patient.

Key Features

- Filter strategically located above the one-way valve to prevent patient secretions from saturating the filter
- Designed for use in EMS, hospital, home care and industrial markets
- Standard mask size accommodates adults and children
- Operation is not affected by extreme temperatures
- Configurations include a hard case for compact storage, oxygen inlet, head strap, gloves and wipes

ORDERING INFORMATION

L62183	Face Mask with one-way valve and filter, 12/case
L64184	Face Mask with one-way valve and filter in hard white case, 12/case
L64184X	Face Mask with one-way valve and filter, gloves & wipes in hard white case, 12/case
L64185	Face Mask with one-way valve and filter, O ₂ inlet & head strap, 12/case
L64186	Face Mask with one-way valve and filter, O ₂ inlet & head strap in hard white case, 12/case
L62050	One-way valve with filter, 50/case



EMERGO EUROPE
Molenstraat 15
2513 BH, The Hague
The Netherlands

L64184
L64186



Flexible CPR Face Shield



L63002 Face Shield
Resuscitator

LSP Flexible CPR Face Shields provide an inexpensive barrier offering protection for the rescuer while administering rescue breathing. This single patient use device incorporates a high-flow one-way valve and a 3M Filtrete™ air filter with 99.8% bacteria filtration efficiency and 99.6% viral filtration efficiency.

Key Features

- Filter strategically located above the one-way valve to prevent patient secretions from saturating the filter
- Directions for use are printed onto the face shield
- Designed for EMS, hospital, home care and industrial markets for both adult and child
- Operation is not affected by extreme temperature

ORDERING INFORMATION

L63002 Flexible CPR Face Shield with one-way valve and filter, 50/case



Allied Healthcare Products, Inc. ▪ 1720 Sublette Avenue ▪ St. Louis, MO 63110 USA
314/771-2400 ▪ 800/444-3954 ▪ FAX: 314/771-0650 ▪ www.alliedhpi.com

U.S.
Toll Free: 800/444-3940
FAX: 800/477-7701

Canada & Mexico
Toll Free: 800/446-0552
FAX: 800/246-6201

International
Telephone: 314/268-1683
FAX: 314/771-5183

ISO 13485: 2003

©Allied Healthcare Products, Inc.®

All specifications are nominal and
subject to change without notice.

Form 90-00-0030
February 2014