

Gomco® Circumcision Clamps

The High Quality Leader in Circumcision Clamps.

Designed to be easy to use and provide a safe and consistent result, these chrome-plated circumcision clamps come in a variety of sizes for different glans diameters from newborn to adult.

Circumcision of the Newborn

With newborns, the clamp generally does not require sutures, rarely produces bleeding, and consistently leaves a clean-cut incision which uneventfully heals in 36 hours (as the mucous membrane and skin become securely united), leaving a clean and natural appearance.

Adult Circumcision

The adult circumcision clamp is a modified version of the infant instrument. Different anesthetics and antiseptics are required for adult circumcision.



02-00-0500 - Extra Small




02-04-0502 - Youth



02-08-0506 - Adult

Gomco® Circumcision Clamps

Part Number	Model	Size	Bottom Flare of Cone Diameter/ Glans Diameter	
02-00-0500	500S	Extra Small	7/16"	1.1 cm
02-01-0500	500	Newborn	1/2"	1.3 cm
02-02-0501	501S	Child	9/16"	1.45 cm
02-03-0501	501	Child	5/8"	1.6 cm
02-04-0502	502	Youth	13/16"	2.1 cm
02-05-0503	503	Adult	1"	2.6 cm
02-06-0504	504	Adult	1-1/8"	2.9 cm
02-07-0505	505	Adult	1-1/4"	3.2 cm
02-08-0506	506	Adult	1-3/8"	3.5 cm



02-08-0506 - Adult

Replacement Parts							
Bell			Plate				
Part Number	Model	Size	Part Number	Model	Size	Part Number	
02-90-0001	500S	Extra Small	02-90-0010	500S	Extra Small	02-90-0019	Arm (used on all models)
02-90-0002	500	Newborn	02-90-0011	500	Newborn	02-90-0020	Nut (used on all models)
02-90-0003	501S	Child	02-90-0012	501S	Child		
02-90-0004	501	Child	02-90-0013	501	Child		
02-90-0005	502	Youth	02-90-0014	502	Youth		
02-90-0006	503	Adult	02-90-0015	503	Adult		
02-90-0007	504	Adult	02-90-0016	504	Adult		
02-90-0008	505	Adult	02-90-0017	505	Adult		
02-90-0009	506	Adult	02-90-0018	506	Adult		

To order Circumcision Clamps or replacement parts, contact your Allied Healthcare sales representative or call Customer Service.



Allied Healthcare Products, Inc. • 1720 Sublette Avenue • St. Louis, MO 63110 USA
314/771-2400 • 800/444-3954 • FAX: 314/771-0650 • www.alliedhpi.com

U.S.
Toll Free: 800/444-3940
FAX: 800/477-7701

Canada & Mexico
Toll Free: 800/446-0552
FAX: 800/246-6201

International
Telephone: 314/268-1683
FAX: 314/771-5183

ISO 13485: 2003

©Allied Healthcare Products, Inc.®

All specifications are nominal and subject to change without notice.

Form 21-00-0003
July 2013